12Biosystems®

Impact of Blood Culture Independent Rapid Diagnostics and Early Identification of Sepsis Causing Pathogens

Aparna Ahuja, MD¹ Devon McCleary, PharmD¹, Elizabeth Olivas¹

¹T2 Biosystems, Lexington, MA, USA

Background

Sepsis contributes between one-third to one-half of all hospital deaths in the USA and accounts for a significant part of the overall costs of healthcare. The incidence of severe sepsis and septic shock in developed countries has been reposted as increasing steadily over recent decades.¹

A study by Pruinelli, L et al., demonstrated that no delay is safe for patients with severe sepsis and septic shock. The longer the delay, the more harm. For example, a delay of 50 minutes, increases the patient's risk of mortality by 0.028 (CI = 0.004, 0.058) suggesting that a delay of 50 minutes

Impact on Lab Results

Figure 3. Culture-Independent Molecular Diagnostic vs Blood Culture Workflow



Impact on Patient Care Results

Table 3. T2Candida: Economic and Clinical Outcomes				
Steuber et al., 2019 ¹⁹	Patch et al., 2018 ²⁰	Kenney et al., 2016 ²¹	Patel et al., 2016 ²²	
 Antifungal optimization occurred in 54% of patients with antifungal orders at time of T2Candida (T2C) test Antifungal therapy was avoided in 60.4% of negative cases Patients with negative T2C results had significantly fewer days of antifungal therapy compared to those with positive results (4.9±6.3 vs 10±10 days, respectively; p=0.03) 	 Faster time to appropriate therapy post-T2C implementation (34 hours to 6 hours, p=0.00147) Unnecessary antifungal therapy avoided in 42.8% and discontinued after 1 dose in another 15.6% of patients with negative T2C results Average length of stay per patient reduced by 8 days Average net antifungal savings of ~\$280 for every patient tested 	 Study demonstrated \$2.3MM in annual hospital savings Reduced median ICU length of stay per patient by 7 days (p=0.009) Reduction in total length of stay by 4 days/patient (p=0.164) 75% of T2C negative patients had antifungals discontinued or deescalated 	 83% of patients who tested positive received appropriate therapy within 6 hours of the blood draw and 100% in under 9 hours Therapy was discontinued for 100% of the patients who tested negative 	

significantly increases the risk of mortality.²

Prompt treatment with targeted antibiotics affects both the financial impact and the clinical outcome of blood stream infection (BSI). Every hour gained in initiating the correct antimicrobial therapy significantly increases the probability of patient survival. Although many infections can be detected after 24 to 48 hours, it may take up to 5 days of incubations to capture most slow-growing bacteria and fungi associated with BSI, and the antibiotic susceptibility determination require an additional 6 to 24 hours. Therefore, a targeted therapy is usually initiated after several days of empirical treatment with broad-spectrum antibiotics. The empiric use of antibiotics results in a 15 to 30% rate of inappropriate treatment, which is associated with a 2 to 5-fold increase in the mortality risk of septic patients and a contributing factors in the recent increases in antimicrobial resistance. ³⁻⁶

T2MR Technology Method

T2 Magnetic Resonance (T2MR) technology from T2Biosystems employs culture-independent testing, providing species identification directly from whole blood samples in 3 to 5 hours without the wait for a positive blood culture.

Faster Time to Species Identification

Table 1: Comparison of T2Bacteria and T2Candida Panels and Blood Culture Time to Result

	Positive Results	T2Candida Panel (hrs.)	Blood Culture (hrs.)	Differenc (hrs.)
Mylonakis et al., 2015 ⁷	Time to Species ID [mean ±SD]	4.4 ± 1	129.9 ± 26.3	103.6
Wilson et al., 2017 ¹³	Time to Species ID [median (IQR)] [‡]	11 (6, 23)	36 (25, 46)	25

	Positive Results	T2Bacteria Panel (hrs.)	Blood Culture (hrs.)	Differen (hrs.)
Nguyen et al., 2019 ⁸	Time to Species ID [mean ± SD]	3.61 ± 0.2 to 7.70 ± 1.38 [†]	71.7 ± 39.3 (with MALDI-TOF or VITEK2)	68.09
De Angelis et al., 2018 ¹⁴	Time to Species ID [mean ± SD]	5.5 ± 1.4	25.2 ± 15.2 (with MALDI-TOF and Film Array BCID)	19.7
†Depending on how many samples were loaded for testing ‡Data presented as median (25 th percentile-75 th percentile)				

Table 4. T2Bacteria: Clinical Outcomes				
Seitz et al., 2019 ²³	Horowitz et al., 2020 ²⁴	Walsh et al., 2019 ²⁵	Voight et al., 2020 ²⁶	
Prospective, single center study of newly admitted patients to Infectious Diseases department (n=44) T2Bacteria (T2B) detected more pathogens than blood culture (BC): 41% (9 of 22) vs 14% (3 of 22) of samples, respectively Faster time to species ID (difference of 59.3 hrs; p 0.01) and faster time to targeted antibiotic therapy (median 6.6 hrs T2B+BC vs 77.7 hrs BC alone) Length of stay was shorter in the T2B+BC group (10.6 days) vs BC alone group (13 days)	 Prospective cohort analysis of Hematopoietic Stem Cell Transplantation patients with febrile neutropenia (n=39, 59 tests) 29 patients had early de- escalation of empiric anti- pseudomonal therapy Total of 124 days of anti- pseudomonal therapy saved; (5.2 days saved per patient) 	 Prospective cohort analysis of Hematologic Malignancy/ HSCT patients (n=94) Faster time to species ID with T2Bacteria 3.7h vs 12.5h with standard blood culture methods (MALDI/BioFire), (p=0.002) T2Bacteria could have potentially influenced care and provided an opportunity to place (T2+/BC-) patients on effective therapy faster than with culture dependent methods 	 Retrospective evaluation of a subset of Emergency Department patients included in the T2Bacteria Pivotal Trial (n=137) Faster time to species ID with T2B (mean 56.6 hours faster) compared to BC T2Bacteria detected 25% more pathogens associated with infection than BC T2Bacteria provided potential opportunity to influence therapy: focus therapy (8 patients), reduce time to species ID (4 patients); reduce time to effective therapy (4 patients) 	

Conclusion



T2Candida and T2Bacteria Panels are the only FDA cleared, commercially available assays for direct-from-blood detection of *Candida* and bacterial species most commonly implicated in BSI that progress to sepsis.

Figure 1. T2DX Workflow



Figure 2. T2Biosystems FDA Cleared Assays

T2Candida®Panel	T2Bacteria® Panel	
Sensitivity: 91% ⁸	Sensitivity: 90%²	
Specificity: 99% ⁸	Specificity: 98%²	
C. albicans	E. faecium	
C. tropicalis	S. aureus	
C. parapsilosis	K. pneumoniae	
C. krusei	P. aeruginosa	
C. glabrata	E. coli	

Improved Diagnosis of Bloodstream Infection

Table 2: T2Candida and T2Bacteria Panels Detect More On-Panel Pathogens Compared to Blood Culture

T2Candida			
Clancy et al., 2019 ¹⁵	Wilson et al., 2017 ¹³		
14-center, prospective trial of hospitalized patients with candidemia (n=152) T2Candida was significantly more likely to be positive that companion blood cultures (cBCs) (45% vs 24%; P<0.0001) Prior antifungal therapy, neutropenia, and <i>C.</i> <i>albicans</i> candidemia were independently associated with T2Candida positivity and T2+/cBC- results (P values <0.05)	 1789 bed, 4 hospital health system; pre-post-test quasi experimental study of hospitalized patients with candidemia (n=161) Candida species were identified in 57% (50/87) patients in the pre-T2Candida group compared to 93% (69/74) in the post-T2Candida group (p<0.001) Median time to appropriate antifungal therapy was reduced from 39 hours to 22 hours (p=0.003) Candida eye involvement was diagnosed in 30% (16/53) in the pre-T2Candida group vs. 12% (4/49) of patients in the post-T2Candida group vs. 12% (P=0.028) 		
T2Ba	cteria		
De Angelis et al., 2018 ¹⁴	Kalligaras at al. 202016		
	Railiger 05 et al., 2020		
1200-bed tertiary-care teaching hospital; prospective observational evaluation of Emergency Department, Infectious Diseases Unit, and ICU patients (n=140samples) T2Bacteria detected 20 positive cases missed by blood culture T2+/BC- results significantly more likely in	 Retrospective evaluation of a subset of patients from the T2Bacteria pivotal trial with discordant T2+/BC- results (n=20 patients, 21 results) 52.5% (11/20) of patients had probable BSI (supported by same organism isolated by at least 1 clinical culture); 19.6% (4/20) possible BSI; 28.5% (6/20) presumptive false-positive 80% (16/20) of patients received active antibiotic 		

SC /			
T2 Biosystems offers an FDA-cleared diagnostic test for the detection of sepsis- causing bacterial and fungal pathogens, directly from whole blood	Species ID within 3-5 hours of first blood draw enables timely therapy when patients need it most	Lab results are often available before the second dose of broad- spectrum antibiotics is delivered	The current blood culture- based standard of care provides species ID results in 2-7 days ¹

Disclosures

AA, DM, and EO, are employees of T2 Biosystems, Inc, the manufacturer of the T2Bacteria and T2Candida Panels.

References

1. Luhr, R., Cao, Y., Söderquist, B. et al. Trends in sepsis mortality over time in randomized sepsis trials: a systematic literature review and meta-analysis of mortality in the control arm, 2002–2016. Crit Care 23, 241 (2019). https://doi.org/10.1186/s13054-019-2528-0 2. Pruinelli L, Westra BL, Yadav P, et al. Delay Within the 3-Hour Surviving Sepsis Campaign Guideline on Mortality for Patients With Severe Sepsis and Septic Shock. Crit Care Med. 2018;46(4):500-505. doi:10.1097/CCM.000000000002949 3. Tassinari M, Zannoli S, Farabegoli P, et al. Rapid diagnosis of bloodstream infections in the critically ill: Evaluation of the broad-range PCR/ESI-MS technology. PLoS One. 2018;13(5):e0197436. Published 2018 May 15. doi:10.1371/journal.pone.0197436 4. Kuti EL, Patel AA, Coleman CI. Impact of inappropriate antibiotic therapy on mortality in patients with ventilator-associated pneumonia and blood stream infection: A meta-analysis. J. Crit. Care. 2008; 23: 91–100. doi: 10.1016/j.jcrc.2007.08.007 5. Kollef MH. Broad-spectrum antimicrobials and the treatment of serious bacterial infections: getting it right up front. Clin Infect Dis. 2008. 47 Suppl 1:S3-6. Ibrahim EH, Sherman G, Ward S, Fraser VJ, Kollef MH. The influence of inadequate antimicrobial treatment of bloodstream infections on patient outcomes in the ICU setting. Chest. 2000; 118(1): 146–155 7. Mylonakis E, Clancy C, Ostrosky-Zeichner L, et al. T2 Magnetic Resonance Assay for the Rapid Diagnosis of Candidemia in Whole Blood: A Clinical Trial. Clinical Infectious Diseases. 2015; 60:892-9 8. Nguyen MH, Clancy CJ, Pasculle AW, et al. Performance of the T2Bacteria Panel for Diagnosing Bloodstream Infections: A Diagnostic Accuracy Study. Ann Intern Med. 2019. Epub 2019/05/15 9. Pfaller M, Kiekema D, Gibbs D, et al. Results from the ARTEMIS DISK Global Antifungal Surveillance Study 1997 to 2007: A 10.5-year analysis of susceptibilities of Candida species to fluconazole and voriconazole as determined by CLSI standardized disk diffusion. J Clin Microbiol 2010b, 48:1366-1377 10. Voight C, Silbert S, Widen R, et al. The T2Bacteria Assay is a sensitive and rapid detector of bacteremia that can be initiated in the emergency department and has potential to favorably influence subsequent therapy. Journal of Emergency Medicine 2020; 58(5): 785-796. 11. Boucher HW, Talbot GH, Bradley JS, et al. Bad bugs, no drugs: no ESKAPE! An update from the Infectious Diseases Society of America. Clin Infect Dis. 2009; 48 (1):1-12. Epub 2008/11/28. 12. Marturano J, Lowery T. ESKAPE pathogens in bloodstream infections are associated with higher cost and mortality but can be predicted using diagnoses upon admission. Open Forum Infectious Diseases 2019; 6(12):1-8. 13. Wilson N, Alangaden G, Tibbetts R, et al. T2 Magnetic Resonance Assay Improves Timely Management of Candidemiia. Journal of Antimicrobial Stewardship 2017; 1(1):12-18. 14. De Angelis G, Posteraro B, De Carolis, et.al. T2Bacteria magnetic resonance assay for the rapid detection of ESKAPEc pathogens directly in whole blood. J Antimicrob Chemother. 2018; 73:iv20-iv26. 15. Clancy C, Pappas P, Vazquez J, etal. Detecting Infections Rapidly and Easily for Candidemia Trial, Part 2 (DIRECT2): A Prospective, Multicenter Study of the T2Candida Panel. Clinical Infectious Diseases 2018.:66: 1678-1686 16. Kalligeros M, Zacharioudakis I, Tansarli G, et al. In-depth analysis of T2Bacteria positive results in patients with concurrent negative blood cultures: a case series. BMC Infectious Diseases 2020; 20(326):1-9. 17. Clancy C, Nguyen M. Finding the "missing 50%" of invasive candidiasis: how nonculture diagnostics will improve understanding of disease spectrum and transform patient care. Clinical infectious diseases 2013; 56(9):1284-1292. 18. Cockerill III, F, Wilson J, Vetter E, et al. Optimal testing parameters for blood cultures. Clinical Infectious Diseases 2004; 38(12):1724-1730. 19. Steuber t, Tucker-Heard G, Edwards J, et al. Utilization and impact of a rapid Candida panel on antifungal stewardship program within a large community hospital. Diagnostic Microbiology and Infectious Disease 2020; 97. 20. Patch M, Weisz E, Cubillos A, et al. Impact of rapid, culture-independent diagnosis of candidaemia and invasive candidiasis in a community health system. J Antimicrob Chemother 2018; 73:4, iv27-iv30. 21. Kenney, R.M., Dwivedi S, Kendal R, et al. Implementation of T2 Magnetic Resonance into the Antimicrobial Stewardship Program Improves Management of Candidemia at HFHS. Poster Presented at: IDWeek, October 2016; New Orleans, LA, USA. 22. Patel F, Young E. Antifungal Prescribing During Initial Implementation of Candidemia Early Detection and Species Identification Testing with T2Candida Panel. Poster Presented at: IDWeek, October 2016; New Orleans LA, USA. 23. Seitz T et al. Evaluation of the clinical impact of the T2MR for the Diagnosis of Blood Stream Infections. Poster presented at: ID Week; October 2019; Washington, DC, USA. 24. Horowitz J et al. Utility of T2Bacteria Panel for Antibiotic Stewardship in Patients Undergoing Hematopoietic Stem Cell Transplantation. Poster presented at: ASTCT; February 2020; Orlando, FL, USA. 25. Walsh T et al. The T2Bacteria Panel is a rapid detector of bacteremia and has potential to guide therapy in patients with hematological malignancies and hematopoietic stem cell transplantation (HSCT): a pilot study of non-culture molecular diagnostics. Poster presented at: ECCMID-EIM; April 2019; Amsterdam, Netherlands, 26. Voight C, Silbert S, Widen R, et al. The T2Bacteria Assay is a sensitive and rapid detector of bacteremia that can be initiated in the emergency department and has potential to favorably influence subsequent therapy. Journal of Emergency Medicine 2020; 58(5): 785-796.



 T2Candida Panel pathogens cover 90% of all Candida blood stream infections⁹

 T2Bacteria Panel pathogens cover 50%-70% of all bacterial blood stream infections^{8,10} and ~90% of ESKAPE pathogens.¹¹

 ESKAPE pathogens are concerning due to their ability to "escape" and survive antibiotic therapy and are associated with increased morbidity and mortality.¹²